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For the Hackuarium Biosafety Level 1 laboratory (P1 Lab), in order to **dispose properly of plates containing genetically modified organisms** (generally bacteria, but also on occasion *C. elegans* strains), they should be inactivated fully, before disposal in the trash.

Follow this procedure for petri dishes:

- 1) Put on personal protective gear: lab coat, gloves and glasses.
- 2) Put plates in a large tray. We have a big white plastic one : 'Clean up container' - generally above the freezer.
- 3) Carefully unseal plates (if sealed) and open them fully.
- 4) Pour 5ml 100% commercial bleach solution* onto each plate. (To use a 15ml conical tube to measure the volume is sufficient.)
- 5) Let sit 5 minutes.
- 6) Agitate a bit, carefully, by hand, still in the tray, checking for dissolving colonies
- 7) Let sit 5 more minutes.
- 8) Seal plates with parafilm and dispose of them in the ordinary trash for incineration.
- 9) Wash out the 'clean up' tray carefully, in back big sink (off the Polymeka space)

For liquid bacterial cultures (not larger than 0.5L):

add 10 ml bleach to 500ml culture (proportionally less for small liquid cultures in tubes), swirl and and let sit 15min before disposing down the drain - wash with fresh water.

***Further notes regarding use of bleach, a toxic CHLORINE COMPOUND**

- Recommended for disinfecting procedures. For disinfection of surfaces, we use ordinary household liquid chlorine bleach, diluted 1:10 with tap water. A 1 to 10 dilution of liquid laundry bleach should contain about 0.5 percent sodium hypochlorite. *For metal surfaces, make sure to wipe down subsequently with 70% ethanol (or alcohol for burners), to avoid corrosion!*
- Bleach decay rates mean the diluted solutions must be made up fresh.
- Bleach is neutralized rapidly in the presence of organic matter.

- Caution should always be taken when using bleach. Use PPE and do not splash it around.

Furthermore:

A 1:10 dilution (5,000 ppm) of bleach has a biocidal effect on *M. tuberculosis*, *S. aureus*, other vegetative bacteria, and HIV after 10-20 minutes. For bacterial spores and mycobacteria, higher concentrations of at least 2,500 ppm (1:5 dilution) are needed.

A 1:50 dilution of chlorine bleach stored at room temperature in a closed plastic container will deteriorate to the equivalent of a 1:100 dilution (500 ppm) after one month, if the container is kept closed (Rutola,W.A., Am. J. Infect Control. 1989;17:1).

[Source](#)